

# Anti-Galectin 3 Rabbit antibody



Catalog Number: 166959

## Product name

Rabbit IgG

Anti-Galectin 3 Rabbit antibody

## Applications

## Specificity

WB, ICC/IF, FC

Human, Rat

## Dilutions

## Antibody description

WB: 1/1000-1/5000

Galectin 3 Rabbit polyclonal antibody

ICC/IF: 1/20-1/50

## Preparation

FC: 1/20

Antigen: A synthetic peptide of human Galectin 3

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type



Western blot detection of Galectin 3 in HeLa, COS7 cell lysates using Galectin 3 Rabbit pAb(1:1000 diluted). Predicted band size: 26KDa. Observed band size: 28KDa.