

# Anti-CD81 Rabbit antibody



Catalog Number: 166963

## Product name

Anti-CD81 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

CD81 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human CD81

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB, IHC-P, ICC/IF, FC


## Dilutions

WB: 1/1000-1/2000

## Validations



Immunohistochemical analysis of paraffin-embedded human tonsil, using CD81 Antibody.



Western blot analysis of CD81 expression in Ramos cell lysate.