

Anti-Fyn Rabbit antibody

Catalog Number: 166995

Product name

Anti-Fyn Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Fyn Rabbit polyclonal antibody

Preparation

Antigen: CTATEPQYQPGENL

Formulation

10mM PBS with 0.02% sodium azide 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot analysis of extracts from Ramos, Jurkat, C6, 3T3 and HeLa cell lysates using Fyn Rabbit pAb (1:1000 diluted). Predicted band size: 59kDa. Observed band size: 45-59kDa.