

Anti-GAB1 Rabbit antibody



Catalog Number: 166998

Product name

Anti-GAB1 Rabbit antibody

Specificity

Human, Rat

Antibody description

GAB1 Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around aa.625~629 (V-E-Y-L-D) derived from Human Gab1.

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB


Dilutions

WB:1/1000

Validations



Western blot analysis of extracts from A431 and MCF7 cells using Gab1 (Ab-627) rabbit pAb (1:1000 diluted). Predicted band size: 110KDa. Observed band size: 110KDa.



Western blot analysis of extracts from A431 and MCF cells using Gab1 (Ab-627) rabbit pAb (1:1000 diluted). Predicted band size: 110KDa. Observed band size: 110KDa.