

Anti-YY1 Rabbit antibody



Catalog Number: 167034

Product name

Anti-YY1 Rabbit antibody

ICC/IF: 1/100-1/500

IP: 1/50-1/100

Specificity

Human, Mouse, Rat

Validations

Antibody description

YY1 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human YY1

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, IP


Dilutions

WB: 1/500-1/3000


IHC: 1/50-1/200



Immunofluorescent analysis of HeLa cells, using YY1 Antibody.



Immunohistochemical analysis of paraffin-embedded human bladder using YY1 Antibody.



Western blot analysis of YY1 expression in (1) HeLa cell lysate; (2) Daudi cell lysate.