

Anti-IL-11RA Rabbit antibody



Catalog Number: 167046

Product name

Anti-IL-11RA Rabbit antibody

WB: 1/1000

IHC: 1/100

Specificity

Human

ICC/IF: 1/50

FC: 1/20-1/100


Antibody description

IL-11RA Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human IL11RA



Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Western blot detection of IL-11RA in K562 cell lysates using IL-11RA Rabbit pAb(1:1000 diluted).Predicted band size:45kDa.Observed band size:45kDa.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions



Immunohistochemistry of IL11RA in paraffin-embedded Human colon cancer tissue using IL11RA Rabbit pAb at dilution 1/50