

Anti-FACL4 Rabbit antibody



Catalog Number: 167232

Product name

Anti-FACL4 Rabbit antibody

IHC: 1/50

ICC/IF: 1/50

Specificity

Human, Mouse, Rat

FC: 1/20

IP: 1/20


Antibody description

FACL4 Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human FACL4



Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of FACL4 in K562,C6,3T3,HeLa cell lysates using FACL4 Rabbit pAb(1:1000 diluted).Predicted band size:79kDa.Observed band size:79kDa.

Clonality


Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP



Dilutions

WB: 1/5000

Immunohistochemistry of FACL4 in paraffin-embedded Human lung cancer tissue using FACL4 Rabbit pAb at dilution 1/50