

Anti-NEK9 Rabbit antibody



Catalog Number: 167238

Product name

Anti-NEK9 Rabbit antibody

WB: 1/2000

IHC: 1/20

Specificity

Human, Mouse, Rat

ICC/IF: 1/50

FC: 1/20


Antibody description

NEK9 Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human NEK9



Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Western blot detection of NEK9 in HeLa cell lysates using NEK9 Rabbit pAb(1:1000 diluted).Predicted band size:107kDa.Observed band size:107kDa.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions



Immunohistochemistry of NEK9 in paraffin-embedded Human colon cancer tissue using NEK9 Rabbit pAb at dilution 1/20