

# Anti-PSS-2 Rabbit antibody



Catalog Number: 167249

## Product name

Anti-PSS-2 Rabbit antibody

## Specificity

Human

## Antibody description

PSS-2 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human PTDSS2

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, FC, IP


## Dilutions

WB: 1/1000-1/2000

FC: 1/20

IP: 1/20

## Validations



Western blot detection of PTDSS2 in K562, HeLa cell lysates using PTDSS2 Rabbit pAb(1:1000 diluted). Predicted band size: 56kDa. Observed band size: 56kDa.