

# Anti-LOXL2 Rabbit antibody



Catalog Number: 167252

## Product name

Anti-LOXL2 Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB, ICC/IF, FC

## Antibody description

LOXL2 Rabbit polyclonal antibody

## Dilutions

WB: 1/500-1/2000

ICC/IF: 1/50-1/200

## Preparation

Antigen: A synthesized peptide derived from human LOXL2

FC: 1/50

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Validations



Western blot analysis of LOXL2 expression in MCF7 cell lysate.