

Anti-LOXL2 Rabbit antibody



Catalog Number: 167252

Product name

Anti-LOXL2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

LOXL2 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human LOXL2

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/500-1/2000

ICC/IF: 1/50-1/200

FC: 1/50

Validations



Western blot analysis of LOXL2 expression in MCF7 cell lysate.