

Anti-LDHA Rabbit antibody



Catalog Number: 167307

Product name

WB: 1/1000

Anti-LDHA Rabbit antibody

IHC: 1/50

Specificity

ICC/IF: 1/50

Human, Mouse, Rat

FC: 1/20

Antibody description

IP: 1/20

LDHA Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Lactate Dehydrogenase

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

Validations



Western blot detection of LDHA in K562, Rat Brain, C6, 3T3, HeLa cell lysates using LDHA Rabbit pAb(1:1000 diluted). Predicted band size: 37kDa. Observed band size: 37kDa.



Immunocytochemistry of Lactate Dehydrogenase(green) in HeLa cells using Lactate Dehydrogenase Rabbit pAb at dilution 1/50, and DAPI(blue)