

# Anti-IL-18 Rabbit antibody



Catalog Number: 167312

## Product name

Anti-IL-18 Rabbit antibody

## Specificity

Human

## Antibody description

IL-18 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human IL18

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC


## Dilutions

WB: 1/1000

IHC: 1/200

FC: 1/20

## Validations



Western blot detection of IL18 in HeLa cell lysates using IL18 Rabbit pAb(1:1000 diluted).Predicted band size:22kDa.Observed band size:17kDa.