

# Anti-Bcl-XL Rabbit antibody



Catalog Number: 167327

## Product name

Anti-Bcl-XL Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Bcl-XL Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Bcl-XL

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC, IP

## Dilutions


WB: 1/1000

IHC: 1/200

FC: 1/20

IP: 1/30

## Validations



Western blot detection of Bcl-xL in HeLa, Jurkat, 3T3, PC12 cell lysates using Bcl-xL Rabbit pAb(1:1000 diluted). Predicted band size: 26KDa. Observed band size: 30KDa.