

Anti-RNA Helicase A Rabbit antibody



Catalog Number: 167348

Product name

Rabbit IgG

Anti-RNA Helicase A Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/2000-1/10000

RNA Helicase A Rabbit polyclonal antibody

ICC/IF: 1/50

Preparation


FC: 1/20

Antigen: A synthetic peptide of human RNA Helicase A

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Western blot detection of RNA Helicase A in C6 cell lysates using RNA Helicase A Rabbit pAb(1:1000 diluted). Predicted band size:141KDa.Observed band size:141KDa.

Ig Type