

Anti-Nrf2 Rabbit antibody



Catalog Number: 167350

Product name

Anti-Nrf2 Rabbit antibody

Specificity

Human

Antibody description

Nrf2 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Nrf2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP, ChIP

Dilutions

WB: 1/1000


ICC/IF: 1/50

FC: 1/20

IP: 1/20

ChIP: 1/20

Validations



Western blot detection of Nrf2 in HepG2 cell lysates using Nrf2 Rabbit pAb(1:1000 diluted).Predicted band size:68KDa.Observed band size:100KDa.