

Anti-Bcl-XL Rabbit antibody



Catalog Number: 167354

Product name

Anti-Bcl-XL Rabbit antibody

IHC: 1/200

ICC/IF: 1/200

Specificity

Human, Mouse, Rat

FC: 1/20

IP: 1/20

Antibody description

Bcl-XL Rabbit polyclonal antibody

Validations

Preparation

Antigen: Recombinant protein of human Bcl-XL

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/5000



Western blot detection of Bcl-XL in Rat Brain, C6, 3T3, HeLa cell lysates using Bcl-XL Rabbit pAb(1:1000 diluted). Predicted band size: 26kDa. Observed band size: 30kDa.



Western blot detection of Bcl-XL in K562, Rat Brain, 3T3 cell lysates using Bcl-XL Rabbit pAb(1:1000 diluted). Predicted band size: 26kDa. Observed band size: 30kDa.