

Anti-Histone H2A.X Rabbit antibody



Catalog Number: 167363

Product name

Anti-Histone H2A.X Rabbit antibody

WB: 1/1000-1/5000

IHC: 1/200

Specificity

Human, Mouse, Rat

ICC/IF: 1/200

FC: 1/20

IP: 1/20

Antibody description

Histone H2A.X Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Histone H2A.X

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

Validations



Western blot detection of Histone H2A.X in 3T3, HeLa cell lysates using Histone H2A.X Rabbit pAb(1:1000 diluted). Predicted band size: 15kDa. Observed band size: 15kDa.



Immunohistochemistry of Histone H2A.X in paraffin-embedded Human colon cancer tissue using Histone H2A.X Rabbit pAb at dilution 1/50