

Anti-Peroxiredoxin 2/PRP Rabbit antibody



Catalog Number: 167399

Product name

Anti-Peroxiredoxin 2/PRP Rabbit antibody

WB: 1/1000-1/5000

IHC: 1/50-1/200

Specificity

Human, Mouse, Rat

ICC/IF: 1/20-1/100

FC: 1/20-1/100

Antibody description

Peroxiredoxin 2/PRP Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human Peroxiredoxin 2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal


Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, FC

Dilutions



Western blot detection of Peroxiredoxin 2/PRP in K562 cell lysates using Peroxiredoxin 2/PRP Rabbit pAb(1:1000 diluted). Predicted band size:22kDa. Observed band size:22kDa.



Immunohistochemistry of Peroxiredoxin 2 in paraffin-embedded Human breast cancer tissue using Peroxiredoxin 2 Rabbit pAb at dilution 1/50