

Anti-Nrf2 Rabbit antibody

Catalog Number: 167417

Product name

Anti-Nrf2 Rabbit antibody

Specificity

Human, Mouse

Antibody description

Nrf2 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human Nrf2

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP, ChIP

Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/200


ICC/IF: 1/50-1/200

IP: 1/50

Validations



Western blot analysis of Nrf2 expression in HepG2



cell lysate.

Immunohistochemical analysis of paraffin-embedded human kidney, using Nrf2 Antibody.



Immunofluorescent analysis of HeLa cells, using Nrf2 Antibody.