

Anti-Hsp27 Rabbit antibody



Catalog Number: 167438

Product name

Anti-Hsp27 Rabbit antibody

WB: 1/1000

IHC: 1/50

Specificity

Human

ICC/IF: 1/20

Antibody description

Hsp27 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Hsp27

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF

Dilutions

Validations



Western blot detection of Hsp27 in HeLa cell lysates using Hsp27 Rabbit pAb(1:1000 diluted).Predicted band size:23kDa.Observed band size:27kDa.



Immunohistochemistry of Hsp27 in paraffin-embedded Human tonsil using Hsp27 Rabbit pAb at dilution 1/50