

Anti-COX2/cyclooxygenase2 Rabbit antibody



Catalog Number: 167457

Product name

Anti-COX2/cyclooxygenase2 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IHC-P, ICC/IF, ELISA

Antibody description

COX2/cyclooxygenase2 Rabbit polyclonal antibody

Dilutions

WB: 1/500 - 1/2000

IHC: 1/100 - 1/300

IF: 1/200 - 1/1000

ELISA: 1/20000

Preparation

Antigen: The antiserum was produced against synthesized peptide derived from human Cox2. AA range:555-604

Validations

Formulation


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal



Western Blot analysis of COx-2 in A549 cells using Cox-2 Polyclonal Antibody at dilution 1/2000