

# Anti-p53 Rabbit antibody



Catalog Number: 167466

## Product name

Anti-p53 Rabbit antibody

## Specificity

Human, Mouse, Rat, Monkey

## Antibody description

p53 Rabbit polyclonal antibody

## Preparation

Antigen: The antiserum was produced against synthesized peptide derived from human p53. AA range:10-59

## Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, ELISA

## Dilutions


WB: 1/500 - 1/2000

IHC: 1/100 - 1/300

IF: 1/200 - 1/1000

ELISA: 1/10000

## Validations



Western Blot analysis of Jurkat cells using p53



Polyclonal Antibody

Western Blot analysis of COS7 cells using p53  
Polyclonal Antibody