

Anti-ATF2 Rabbit antibody



Catalog Number: 167468

Product name

Anti-ATF2 Rabbit antibody

Specificity

Human

Antibody description

ATF2 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human ATF2

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/100


ICC/IF: 1/50-1/100

IP: 1/30

Validations



Western blot analysis of ATF2 expression in HeLa cell lysate.



Immunofluorescent analysis of HeLa cells, using ATF2 Antibody .