

Anti-p21 Rabbit antibody

Catalog Number: 167522

Product name

Anti-p21 Rabbit antibody

ICC/IF: 1/50-1/200

IP: 1/50


Specificity

Human

Validations

Antibody description

p21 Rabbit polyclonal antibody




Preparation

Antigen: A synthesized peptide derived from human p21

Immunofluorescent analysis of MCF7 cells, using p21 Antibody .

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Immunohistochemical analysis of paraffin-embedded human colon, using p21 Antibody.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, IP

Dilutions

WB: 1/500-1/2000

IHC: 1/50-1/200



Western blot analysis of p21 in (1) MCF-7 cell lysate; (2) LnCaP cell lysate.