

Anti-GC1q R Rabbit antibody



Catalog Number: 167540

General Information

Product name:

Anti-GC1q R Rabbit antibody

Immunogen:

Applications:

WB, IP

Specificity:

Human, Mouse, Rat

Antibody description:

GC1q R Rabbit polyclonal antibody

Properties

Preparation:

Antigen: A synthetic peptide of human C1QBP

Formulation:

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage:

Store at -20°C. Stable for 12 months from date of receipt.

Clonality:

Polyclonal


Ig Type:Rabbit IgG

Applications

WB: 1/1000-1/5000

IP: 1/20

Validations



Western blot detection of GC1q R in K562,Rat Brain,C6,3T3,Hela cell lysates using GC1q R Rabbit pAb(1:1000 diluted).Predicted band size:31kDa.Observed band size:31kDa.