

# Anti-Smad3 Rabbit antibody



Catalog Number: 167551

## Product name

Anti-Smad3 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Smad3 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human Smad3

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions


WB: 1/1000

IP: 1/20

## Validations



Western blot detection of Smad3 in K562,CHO-K1 cell lysates using Smad3 Rabbit pAb(1:1000 diluted).Predicted band size:48KDa.Observed



band size:52KDa.

Western blot detection of Smad3 in Jurkat cell lysates using Smad3 Rabbit pAb(1:1000 diluted).Predicted band size:48KDa.Observed band size:52KDa.