

# Anti-Amyloid Precursor Protein Rabbit antibody



Catalog Number: 167552

## Product name

Anti-Amyloid Precursor Protein Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IP

## Antibody description

Amyloid Precursor Protein Rabbit polyclonal antibody

## Dilutions


WB: 1/1000

IP: 1/20

## Preparation

Antigen: A synthetic peptide of human Amyloid Precursor Protein

## Validations



## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of Amyloid Precursor Protein in HeLa cell lysates using Amyloid Precursor Protein Rabbit pAb(1:500 diluted). Predicted band size: 87kDa. Observed band size: 100kDa.

## Clonality

Polyclonal