

Anti-Amyloid Precursor Protein Rabbit antibody



Catalog Number: 167552

Product name

Anti-Amyloid Precursor Protein Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB, IP

Antibody description

Amyloid Precursor Protein Rabbit polyclonal antibody

Dilutions

WB: 1/1000

IP: 1/20

Preparation

Antigen: A synthetic peptide of human Amyloid Precursor Protein

Validations

Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal



Western blot detection of Amyloid Precursor Protein in HeLa cell lysates using Amyloid Precursor Protein Rabbit pAb(1:500 diluted). Predicted band size: 87kDa. Observed band size: 100kDa.