

Anti-IL-1RAcP Rabbit antibody



Catalog Number: 167571

Product name

Anti-IL-1RAcP Rabbit antibody

Specificity

Human

Antibody description

IL-1RAcP Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human IL-1RAcP

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, FC, IP


Dilutions

WB: 1/1000

FC: 1/20

IP: 1/20

Validations



Western blot detection of IL-1RAcP in HL-60,U251,U87-MG cell lysates using IL-1RAcP Rabbit mAb(1:1000 diluted).Predicted band size:65kDa.Observed band size:78kDa.