

# Anti-NEDD4 Rabbit antibody



Catalog Number: 167594

## Product name

Anti-NEDD4 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

NEDD4 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human NEDD4

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB, FC, IP


## Dilutions

WB: 1/2000-1/10000

FC: 1/20

IP: 1/20

## Validations



Western blot detection of NEDD4 in HeLa,A549,HL-60,U2OS,U87-MG cell lysates using NEDD4 Rabbit mAb(1:1000 diluted).Predicted band size:149kDa.Observed band size:115kDa.