

# Anti-Glucocorticoid Receptor Rabbit antibody



Catalog Number: 167618

## Product name

Anti-Glucocorticoid Receptor Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Glucocorticoid Receptor Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Glucocorticoid Receptor

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/2000-1/10000

FC: 1/20-1/100

## Validations



Western blot detection of Glucocorticoid Receptor in CHO-K1,C6,Rat Brain,Ramos cell lysates using Glucocorticoid Receptor Rabbit mAb(1:1000 diluted).Predicted band size:86kDa.Observed band size:86kDa.