

# Anti-OR1A1 Rabbit antibody



Catalog Number: 167621

## Product name

Anti-OR1A1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

OR1A1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human OR1A1

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications


WB, IP

## Dilutions

WB: 1/1000-1/2000

IP: 1/20

## Validations



Western blot detection of OR1A1 in CHO-K1, C6, Jurkat cell lysates using OR1A1 Rabbit mAb(1:1000 diluted). Predicted band size: 35kDa. Observed band size: 35kDa.