

Anti-phospho-Cyclin B1 (Ser126) Rabbit antibody



Catalog Number: 167625

Product name	Monoclonal
Anti-phospho-Cyclin B1 (Ser126) Rabbit antibody	Ig Type
Specificity	Rabbit IgG
Human	Applications
Antibody description	WB, IP
phospho-Cyclin B1 (Ser126) Rabbit monoclonal antibody	Dilutions
Preparation	WB: 1/1000
	IP: 1/20

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser126 of human Cyclin B1

Formulation


50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Validations



Western blot detection of Cyclin B1 (phospho-Ser126) in HeLa,Ramos using Cyclin B1 (phospho-Ser126) Rabbit mAb(1:1000 diluted)