

Anti-LIG1 Rabbit antibody



Catalog Number: 167660

Product name

Anti-LIG1 Rabbit antibody

Specificity

Human, Rat

Antibody description

LIG1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human LIG1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000

IHC: 1/20

ICC/IF: 1/50

Validations



Western blot detection of LIG1 in HeLa,A549,HL-60 using LIG1 Rabbit mAb(1:1000 diluted)