

# Anti-Clathrin Heavy Chain Rabbit antibody



Catalog Number: 167662

## Product name

Anti-Clathrin Heavy Chain Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Clathrin Heavy Chain Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Clathrin heavy chain

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications


WB, IHC-P

## Dilutions

WB: 1/2000-1/10000

IHC: 1/20

## Validations



Western blot detection of Clathrin heavy chain in HeLa,A549,HL-60,U251,U87-MG cell lysates using Clathrin heavy chain Rabbit mAb(1:1000 diluted).Predicted band size:192kDa.Observed band size:192kDa.