

Anti-CTBP2 Rabbit antibody



Catalog Number: 167693

Product name

Anti-CTBP2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

CTBP2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human CTBP2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/2000-1/10000

IHC: 1/50

ICC/IF: 1/50

Validations



Western blot detection of CTBP2 in HeLa, A549, U2OS, C6 cell lysates using CTBP2 Rabbit mAb(1:1000 diluted). Predicted band size: 49kDa. Observed band size: 49kDa.