

Anti-Cofilin Rabbit antibody



Catalog Number: 167694

Product name

Anti-Cofilin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Cofilin Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Cofilin

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

FC: 1/200

Validations



Western blot detection of Cofilin in A549, MCF-7, 293T cell lysates using Cofilin Rabbit mAb(1:1000 diluted). Predicted band size: 19kDa. Observed band size: 19kDa.