

# Anti-Hex Rabbit antibody



Catalog Number: 167700

## Product name

Rabbit IgG

Anti-Hex Rabbit antibody

## Applications

## Specificity

WB, IP

Human

## Dilutions

## Antibody description

WB: 1/1000

Hex Rabbit monoclonal antibody

IP: 1/20

## Preparation

Antigen: A synthetic peptide of human Hex

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage


Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

## Validations



Western blot detection of Hex in Jurkat,HeLa cell lysates using Hex Rabbit mAb(1:1000 diluted).Predicted band size:30kDa.Observed band size:30kDa.