

Anti-phospho-TAOK1(Ser181)/TAOK2 (Ser181)/TAOK3 (Ser177) Rabbit antibody



Catalog Number: 167723

Product name

Anti-phospho-TAOK1(Ser181)/TAOK2 (Ser181)/TAOK3 (Ser177) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

phospho-TAOK1(Ser181)/TAOK2 (Ser181)/TAOK3 (Ser177) Rabbit monoclonal antibody

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser181 of human TAOK1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/50

Validations



Western blot detection of TAOK1 (phospho-Ser181)/TAOK2 (phospho-Ser181)/TAOK3 (phospho-Ser177) in HeLa, A549, HL-60, U2OS, C6 cell lysates using TAOK1 (phospho-Ser181)/TAOK2 (phospho-Ser181)/TAOK3 (phospho-Ser177) Rabbit mAb(1:1000 diluted). Predicted band size: 105, 116, 138 kDa. Observed band size: 105, 116, 138 kDa.