

# Anti-Poliovirus Receptor/PVR Rabbit antibody



Catalog Number: 167761

## Product name

Anti-Poliovirus Receptor/PVR Rabbit antibody

## Specificity

Human

## Antibody description

Poliovirus Receptor/PVR Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Poliovirus Receptor/PVR

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB, IHC-P, FC


## Dilutions

WB: 1/1000

IHC: 1/50

FC: 1/20

## Validations



Western blot detection of Poliovirus Receptor/PVR in HeLa, A549, U2OS cell lysates using Poliovirus Receptor/PVR Rabbit mAb(1:1000 diluted). Predicted band size: 45kDa. Observed band size: 60-80kDa.