

# Anti-phospho-p53 (Ser6) Rabbit antibody



Catalog Number: 167771

## Product name

Rabbit IgG

Anti-phospho-p53 (Ser6) Rabbit antibody

## Applications

## Specificity

WB, IHC-P, IP

Human

## Dilutions

## Antibody description

WB: 1/1000

phospho-p53 (Ser6) Rabbit monoclonal antibody

IHC: 1/20

## Preparation

IP: 1/20

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser6 of human p53

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of phospho-p53 (Ser6) in A549, MCF-7 cell lysates using phospho-p53 (Ser6) Rabbit mAb (1:1000 diluted). Predicted band size: 44 kDa. Observed band size: 53 kDa.