

Anti-TPMT Rabbit antibody



Catalog Number: 167866

Product name

Anti-TPMT Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TPMT Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human TPMT

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/50

FC: 1/20-1/100

IP: 1/20

Validations



Western blot detection of TPMT in HeLa,A549,HL-60 cell lysates using TPMT Rabbit mAb(1:1000 diluted).Predicted band size:28kDa.Observed band size:28kDa.