

Anti-Poliovirus Receptor/PVR Rabbit antibody



Catalog Number: 167886

Product name

Anti-Poliovirus Receptor/PVR Rabbit antibody

Specificity

Human

Antibody description

Poliovirus Receptor/PVR Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Poliovirus Receptor/PVR

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/50

Validations



Western blot detection of Poliovirus Receptor/PVR in HeLa,A549,U2OS cell lysates using Poliovirus Receptor/PVR Rabbit mAb(1:1000 diluted).Predicted band size:45kDa.Observed band size:60-80kDa.