

# Anti-Rad51 Rabbit antibody



Catalog Number: 167906

## Product name

Anti-Rad51 Rabbit antibody

## Specificity

Human

## Antibody description

Rad51 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Rad51

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

IHC-P, ICC/IF, FC


## Dilutions

IHC: 1/50

ICC/IF: 1/20

FC: 1/50

## Validations



Western blot detection of Rad51 in A549,HL-60,U2OS cell lysates using Rad51 Rabbit mAb(1:1000 diluted).Predicted band size:37kDa.Observed band size:37kDa.