

Anti-Tuberin Rabbit antibody



Catalog Number: 167911

Product name

Anti-Tuberin Rabbit antibody

Specificity

Human, Mouse

Antibody description

Tuberin Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human Tuberin

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of Tuberin in HeLa,A549,HL-60,U2OS cell lysates using Tuberin Rabbit mAb(1:1000 diluted).Predicted band size:201kDa.Observed band size:201kDa.