

Anti-Eph receptor B1/NET Rabbit antibody



Catalog Number: 167913

Product name

Anti-Eph receptor B1/NET Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IP

Antibody description

Eph receptor B1/NET Rabbit monoclonal antibody

Dilutions


WB: 1/1000

IP: 1/20

Preparation

Antigen: A synthetic peptide of human Eph receptor B1/NET

Validations



Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Western blot detection of Eph receptor B1/NET in A549,MCF-7,U251 cell lysates using Eph receptor B1/NET Rabbit mAb(1:1000 diluted).Predicted band size:108kDa.Observed band size:135kDa.