

# Anti-Eph receptor B1/NET Rabbit antibody



Catalog Number: 167913

## Product name

Anti-Eph receptor B1/NET Rabbit antibody

## Ig Type

IgG

## Specificity

Human, Mouse, Rat

## Applications

WB, IP

## Antibody description

Eph receptor B1/NET Rabbit monoclonal antibody

## Dilutions

WB: 1/1000

IP: 1/20

## Preparation

Antigen: A synthetic peptide of human Eph receptor B1/NET

## Validations

## Formulation


50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of Eph receptor B1/NET in A549, MCF-7, U251 cell lysates using Eph receptor B1/NET Rabbit mAb(1:1000 diluted). Predicted band size: 108kDa. Observed band size: 135kDa.