

Anti-SNAP25 Rabbit antibody



Catalog Number: 167952

Product name

Anti-SNAP25 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SNAP25 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human SNAP25

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/50

FC: 1/50

IP: 1/20

Validations



Western blot detection of SNAP25 in HL-60,C6 using SNAP25 Rabbit mAb(1:1000 diluted)