

Anti-Cdk9 Rabbit antibody



Catalog Number: 167955

Product name

Anti-Cdk9 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Cdk9 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Cdk9

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/20

FC: 1/50

IP: 1/20

Validations



Western blot detection of Cdk9 in HeLa, MCF-7, U2OS, U251 using Cdk9 Rabbit mAb(1:1000 diluted)